


VORTEX MIXER (NB-105V)



Optional Accessory

| | |
|-----------|------------------------|
| NB-105V01 | Mixing Cup |
| NB-105V02 | Platform Head |
| NB-105V03 | Micro plate Kit |
| NB-105V04 | Multi sample plate Kit |



Features

- Automatic start when vortex is touched on the mixing cup.
- Constant on after starting
- Adjustable RPM up to 3000rpm
- Wide option for accessories
- High quality design for safety
- CE marked

Specification

| Items | Unit | NB-105V |
|----------------|------|----------------------------|
| Speed range | rpm | 0 to 3000rpm |
| Operation Mode | | Touch On & Constant On |
| Material | | Aluminum Die-casting & ABS |
| Head | mm | Rubber Mixing Cup |
| Dimensions | kg | 130(W)x155(D)x155(H)mm 3Kg |
| Power | V/Hz | AC220V, 0.14A, 60Hz |