

Crimping Tools and Decapping tools



- Coating out of the handle, ensure the tools can prevent hull corrosion;
- The head of the tools made from Special alloy materials ensure the long life;
- The handle part use frosted design ensure easy clench;
- Can be adjust the screw nut on the handle to ensure the best sealing performance;
- Crimping tools can use hexagon spanner to adjust the crimp neck height;
- Decapper easy to use, avoid the samples spilling and the glass parts broken;

Crimping tools

Description	Packaging	Cat. No.
Crimper for 8mm Aluminum caps	1 pc. Per box	3.009300.0800
Crimper for 11mm Aluminum caps	1 pc. Per box	3.009300.1100
Crimper for 13mm Aluminum caps	1 pc. Per box	3.009300.1300
Crimper for 20mm Aluminum caps	1 pc. Per box	3.009300.2000

Decapping tools

Description	Packaging	Cat. No.
Decapper for 8mm Aluminum caps	1 pc. Per box	3.009320.0800
Decapper for 11mm Aluminum caps	1 pc. Per box	3.009320.1100
Decapper for 13mm Aluminum caps	1 pc. Per box	3.009320.1300
Decapper for 20mm Aluminum caps	1 pc. Per box	3.009320.2000

Standard screw thread neck vials

These large volume vials is best choice for laboratory to storage the samples.

15-425 standard screw thread neck 8mL vials



15-425 screw thread neck vials

Description	Packaging	Cat. No.
15-425 Thread screw neck vial, 8ml, 61x16.6mm, clear glass, Borosilicate Type I Class A	100pcs. per Carton, 8 Cartons per Case	3.038015.0000
15-425 Thread screw neck vial, 8ml, 61x16.6mm, amber glass, Borosilicate Type I Class B	100pcs. per Carton, 8 Cartons per Case	3.038015.000A
15-425 Thread screw neck vial, 8ml, 61x16.6mm, clear glass, white marking spot and CNW LOGO, Borosilicate Type I Class A	100pcs. per Carton, 8 Cartons per Case	3.038015.00E0
15-425 Thread screw neck vial, 8ml, 61x16.6mm, amber glass, white marking spot and CNW LOGO, Borosilicate Type I Class B	100pcs. per Carton, 8 Cartons per Case	3.038015.00EA

15-425 screw thread neck preassembled caps

Description	Packaging	Cat. No.
Preassembled cap and septa for 15-425 thread screw, PP cap, black, centre hole, Clear silicone/White PTFE, 0.060" thick	100pcs. per PE bag	3.005350.1500
Preassembled cap and septa for 15-425 thread screw, PP cap, black, closed, Clear silicone/White PTFE, 0.060" thick	100pcs. per PE bag	3.005360.1500
Preassembled cap and septa for 15-425 thread screw, PP cap, black, centre hole, Red Butyl/Grey PTFE, 0.060" thick	100pcs. per PE bag	3.0F5350.1500
Preassembled cap and septa for 15-425 thread screw, PP cap, black, closed, Red Butyl/Grey PTFE, 0.060" thick	100pcs. per PE bag	3.0F5360.1500
Preassembled cap and septa for 15-425 thread screw, PP cap, white, closed, F217 line	100pcs. per PE bag	3.W05360.1500

15-425 screw thread neck caps

Description	Packaging	Cat. No.
15-425 Screw thread cap, made from PP, black, centre hole	100pcs. per PE bag	3.005310.1500
15-425 Screw thread cap, made from PP, black, closed	100pcs. per PE bag	3.005320.1500



15-425 screw thread neck septas

Description	Packaging	Cat. No.
Septa for 15-425 screw thread cap, white PTFE only, 0.010" thick	100pcs. per PE bag	3.601010.1500
Septa for 15-425 screw thread cap, Red Butyl/Grey PTFE, 0.060" thick	100pcs. per PE bag	3.606040.150F
Septa for 15-425 screw thread cap, Clear silicone/White PTFE, 0.060" thick	100pcs. per PE bag	3.606050.150G
Septa for 15-425 screw thread cap, Clear silicone/White PTFE, 0.060" thick, ultra low bleed	100pcs. per PE bag	3.606050.15LB




18-400 standard screw thread neck 15mL vials

18-400 screw thread neck vials

Description	Packaging	Cat. No.
18-400 Thread screw neck vial, 15ml, 71x20.6mm, clear glass, Borosilicate Type I Class A	100pcs. per Carton, 10 Cartons per Case	3.316018.0000
18-400 Thread screw neck vial, 15ml, 71x20.6mm, amber glass, Borosilicate Type I Class B	100pcs. per Carton, 10 Cartons per Case	3.316018.000A
18-400 Thread screw neck vial, 15ml, 71x20.6mm, clear glass, white marking spot and CNW LOGO, Borosilicate Type I Class A	100pcs. per Carton, 10 Cartons per Case	3.316018.00E0
18-400 Thread screw neck vial, 15ml, 71x20.6mm, amber glass, white marking spot and CNW LOGO, Borosilicate Type I Class B	100pcs. per Carton, 10 Cartons per Case	3.316018.00EA

18-400 screw thread neck preassembled caps

Description	Packaging	Cat. No.
Preassembled cap and septa for 18-400 thread screw, PP cap, black, centre hole, Clear silicone/White PTFE, 0.060" thick	100pcs. per PE bag	3.005350.1800
Preassembled cap and septa for 18-400 thread screw, PP cap, black, closed, Clear silicone/White PTFE, 0.060" thick	100pcs. per PE bag	3.005360.1800
Preassembled cap and septa for 18-400 thread screw, PP cap, black, centre hole, Red Butyl/Grey PTFE, 0.060" thick	100pcs. per PE bag	3.0F5350.1800
Preassembled cap and septa for 18-400 thread screw, PP cap, black, closed, Red Butyl/Grey PTFE, 0.050" thick	100pcs. per PE bag	3.0F5360.1800




18-400 screw thread neck caps

Description	Packaging	Cat. No.
18-400 Screw thread cap, made from PP, black, centre hole	100pcs. per PE bag	3.005310.1800
18-400 Screw thread cap, made from PP, black, closed	100pcs. per PE bag	3.005320.1800

18-400 screw thread neck septas

Description	Packaging	Cat. No.
Septa for 18-400screw thread cap, white PTFE only, 0.010" thick	100pcs. per PE bag	3.601010.1800
Septa for 18-400 screw thread cap, Red Butyl/Grey PTFE, 0.060" thick	100pcs. per PE bag	3.606040.180F
Septa for 18-400 screw thread cap, Clear silicone/White PTFE, 0.060" thick, ultra low bleed	100pcs. per PE bag	3.606050.18LB
Septa for 18-400 screw thread cap, Clear silicone/White PTFE, 0.100" thick	100pcs. per PE bag	3.610050.180L



24-400 standard screw thread neck 20mL, 40mL, 60mL EPA/VOA vials

- You can use them on below brand instrument: Agilent, Dionex, Shimadzu, Tekmar, Thermo Scientific, Varian.
- Broad range of EPA vials in clear and amber glass
- Volumes of 20mL, 20mL, 40mL and 60mL available

24-400 screw thread neck vials

Description	Packaging	Cat. No.
24-400 Thread screw neck EPA vial, 20ml, 57x27.5mm, clear glass, Borosilicate Type I Class A	100pcs. per Carton, 8 Cartons per Case	3.320024.0000
24-400 Thread screw neck EPA vial, 20ml, 57x27.5mm, amber glass, Borosilicate Type I Class B	100pcs. per Carton, 8 Cartons per Case	3.320024.000A
24-400 Thread screw neck EPA vial, 20ml, 57x27.5mm, clear glass, white marking spot and CNW LOGO, Borosilicate Type I Class A	100pcs. per Carton, 8 Cartons per Case	3.320024.00E0

Standard screw thread neck vials

Vials



Description	Packaging	Cat. No.
24-400 Thread screw neck EPA vial, 20ml, 57x27.5mm, amber glass, white marking spot and CNW LOGO, Borosilicate Type I Class B	100pcs. per Carton, 8 Cartons per Case	3.320024.00EA
24-400 Thread screw neck EPA vial, 20ml, 57x27.5mm, clear glass, Borosilicate type 70	100pcs. per Carton, 8 Cartons per Case	3.320024.7000
24-400 Thread screw neck EPA vial, 20ml, 57x27.5mm, amber glass, Borosilicate type 70	100pcs. per Carton, 8 Cartons per Case	3.320024.700A
24-400 Thread screw neck EPA vial, 30ml, 72.5x27.5mm, clear glass, Borosilicate Type I Class A	100pcs. per Carton	3.330024.0000
24-400 Thread screw neck EPA vial, 30ml, 72.5x27.5mm, amber glass, Borosilicate Type I Class B	100pcs. per Carton	3.330024.000A
24-400 Thread screw neck EPA vial, 30ml, 72.5x27.5mm, clear glass, white marking spot and CNW LOGO, Borosilicate Type I Class A	100pcs. per Carton	3.330024.00E0
24-400 Thread screw neck EPA vial, 30ml, 72.5x27.5mm, amber glass, white marking spot and CNW LOGO, Borosilicate Type I Class B	100pcs. per Carton	3.330024.00EA
24-400 Thread screw neck EPA vial, 30ml, 72.5x27.5mm, clear glass, Borosilicate type 70	100pcs. per Carton	3.330024.7000
24-400 Thread screw neck EPA vial, 30ml, 72.5x27.5mm, amber glass, Borosilicate type 70	100pcs. per Carton	3.330024.700A
24-400 Thread screw neck EPA vial, 40ml, 95x27.5mm, clear glass, Borosilicate Type I Class A	100pcs. per Carton, 4 Cartons per Case	3.340024.0000
24-400 Thread screw neck EPA vial, 40ml, 95x27.5mm, amber glass, Borosilicate Type I Class B	100pcs. per Carton, 4 Cartons per Case	3.340024.000A
24-400 Thread screw neck EPA vial, 40ml, 95x27.5mm, clear glass, white marking spot and CNW LOGO, Borosilicate Type I Class A	100pcs. per Carton, 4 Cartons per Case	3.340024.00E0
24-400 Thread screw neck EPA vial, 40ml, 95x27.5mm, amber glass, white marking spot and CNW LOGO, Borosilicate Type I Class B	100pcs. per Carton, 4 Cartons per Case	3.340024.00EA
24-400 Thread screw neck EPA vial, 40ml, 95x27.5mm, clear glass, Borosilicate type 70	100pcs. per Carton, 4 Cartons per Case	3.340024.7000
24-400 Thread screw neck EPA vial, 40ml, 95x27.5mm, amber glass, Borosilicate type 70	100pcs. per Carton, 4 Cartons per Case	3.340024.700A
24-400 Thread screw neck EPA vial, 60ml, 140x27.5mm, clear glass, Borosilicate Type I Class A	100pcs. per Carton, 4 Cartons per Case	3.360024.0000
24-400 Thread screw neck EPA vial, 60ml, 140x27.5mm, amber glass, Borosilicate Type I Class B	100pcs. per Carton, 4 Cartons per Case	3.360024.000A
24-400 Thread screw neck EPA vial, 60ml, 140x27.5mm, clear glass, white marking spot and CNW LOGO, Borosilicate Type I Class A	100pcs. per Carton, 4 Cartons per Case	3.360024.00E0
24-400 Thread screw neck EPA vial, 60ml, 140x27.5mm, amber glass, white marking spot and CNW LOGO, Borosilicate Type I Class B	100pcs. per Carton, 4 Cartons per Case	3.360024.00EA
24-400 Thread screw neck EPA vial, 60ml, 140x27.5mm, clear glass, Borosilicate type 70	100pcs. per Carton, 4 Cartons per Case	3.360024.7000
24-400 Thread screw neck EPA vial, 60ml, 140x27.5mm, amber glass, Borosilicate type 70	100pcs. per Carton, 4 Cartons per Case	3.360024.700A

24-400 screw thread neck preassembled caps


Description	Packaging	Cat. No.
Preassembled cap and septa for 24-400 thread screw, PP cap, black, centre hole, Clear silicone/White PTFE, 0.100" thick	100pcs. per PE bag	3.005350.2400
Preassembled cap and septa for 24-400 thread screw, PP cap, black, closed, Clear silicone/White PTFE, 0.100" thick	100pcs. per PE bag	3.005360.2400
Preassembled cap and septa for 24-400 thread screw, PP cap, black, centre hole, Red Butyl/Grey PTFE, 0.060" thick	100pcs. per PE bag	3.0F5350.2400
Preassembled cap and septa for 24-400 thread screw, PP cap, black, closed, Red Butyl/Grey PTFE, 0.060" thick	100pcs. per PE bag	3.0F5360.2400
Preassembled cap and septa for 24-400 thread screw, PP cap, white, centre hole, Clear silicone/White PTFE, 0.100" thick	100pcs. per PE bag	3.W05350.2400
Preassembled cap and septa for 24-400 thread screw, PP cap, white, closed, Clear silicone/White PTFE, 0.100" thick	100pcs. per PE bag	3.W05360.2400

Description	Packaging	Cat. No.
Preassembled cap and septa for 24-400 thread screw, PP cap, white, centre hole, Red Butyl/Grey PTFE, 0.100" thick	100pcs. per PE bag	3.WF5350.2400
Preassembled cap and septa for 24-400 thread screw, PP cap, white, closed, Red Butyl/Grey PTFE, 0.060" thick	100pcs. per PE bag	3.WF5360.2400



24-400 screw thread neck caps

Description	Packaging	Cat. No.
24-400 Screw thread cap, made from PP, black, centre hole	100pcs. per PE bag	3.005310.2400
24-400 Screw thread cap, made from PP, black, closed	100pcs. per PE bag	3.005320.2400
24-400 Screw thread cap, made from PP, white, centre hole	100pcs. per PE bag	3.W05310.2400
24-400 Screw thread cap, made from PP, white, closed	100pcs. per PE bag	3.W05320.2400



24-400 screw thread neck septum

Description	Packaging	Cat. No.
Septa for 24-400 screw thread cap, white PTFE only, 0.010" thick	100pcs. per PE bag	3.601010.2400
Septa for 24-400 screw thread cap, Red Butyl/Grey PTFE, 0.060" thick	100pcs. per PE bag	3.606040.240F
Septa for 24-400 screw thread cap, Clear silicone/White PTFE, 0.100" thick	100pcs. per PE bag	3.610050.240L
Septa for 24-400 screw thread cap, Clear silicone/White PTFE, 0.100" thick, ultra low bleed	100pcs. per PE bag	3.610050.24LB

Sampling Bottles

- These large volume sample bottles are so related to environmental samples. According to customer's requirements, we provide the corresponding product.
- Large open mouth bottles are designed for storing the soil, sediments and sludge samples.
- Narrow open mouth bottles are suitable for storing liquid samples.
- All the sampling bottles with cap and the inner liner.
- Brown sample bottles are suitable for storing light-sensitive samples.



Straight bottle

Suitable for storing the soil, sediments and sludge samples.

Description	Packaging	Cat. No.
30mL flint glass jar, screw finish 38/R3, 38/R3 mm black screw cap, material phenolic moulding poeder, mounted waxed KKPVDC-liner	48pcs. Per Carton	3.P14801.0300
60mL flint glass jar, screw finish 51/R3, 51/R3 mm black screw cap, material phenolic moulding poeder, mounted waxed KKPVDC-liner	48pcs. Per Carton	3.P14801.0600
120mL flint glass jar, screw finish 58/R3, 58/R3 mm black screw cap, material phenolic moulding poeder, mounted waxed KKPVDC-liner	12pcs. Per Carton	3.P14801.1200

large open mouth bottle

Suitable for storing the soil, sediments and sludge samples.



Description	Packaging	Cat. No.
Weithalsglaeser, Material:Glas Hydrolytische Klasse III, weiB, 30mL, DIN 32, Verschluesse fuer braune und weiBe Weithalsglaeser, schwarz, PP, LKD, DIN32	48pcs. Per Carton	3.P14802.0300
Weithalsglaeser, Material:Glas Hydrolytische Klasse III, braun, 30mL, DIN 32, Verschluesse fuer braune und weiBe Weithalsglaeser, schwarz, PP, LKD, DIN32	48pcs. Per Carton	3.P14802.030A
Weithalsglaeser, Material:Glas Hydrolytische Klasse III, weiB, 50mL, DIN 32, Verschluesse fuer braune und weiBe Weithalsglaeser, schwarz, PP, LKD, DIN32	24pcs. Per Carton	3.P14802.0500
Weithalsglaeser, Material:Glas Hydrolytische Klasse III, braun, 50mL, DIN 32, Verschluesse fuer braune und weiBe Weithalsglaeser, schwarz, PP, LKD, DIN32	48pcs. Per Carton	3.P14802.050A
Weithalsglaeser, Material:Glas Hydrolytische Klasse III, weiB, 125mL, DIN 40, Verschluesse fuer braune und weiBe Weithalsglaeser, schwarz, PP, LKD, DIN40	12pcs. Per Carton	3.P14802.1250



Description	Packaging	Cat. No.
Weithalsglaeser, Material:Glas Hydrolytische Klasse III, braun, 125mL, DIN 40, Verschluesse fuer braune und weiBe Weithalsglaeser, schwarz, PP, LKD, DIN40	12pcs. Per Carton	3.P14802.125A
Weithalsglaeser, Material:Glas Hydrolytische Klasse III, weiB, 250mL, DIN 55, Verschluesse fuer braune und weiBe Weithalsglaeser, schwarz, PP, LKD, DIN55	12pcs. Per Carton	3.P14802.2500
Weithalsglaeser, Material:Glas Hydrolytische Klasse III, braun, 250mL, DIN 55, Verschluesse fuer braune und weiBe Weithalsglaeser, schwarz, PP, LKD, DIN55	12pcs. Per Carton	3.P14802.250A
Weithalsglaeser, Material:Glas Hydrolytische Klasse III, weiB, 500mL, DIN 55, Verschluesse fuer braune und weiBe Weithalsglaeser, schwarz, PP, LKD, DIN55	6pcs. Per Carton	3.P14802.5000
Weithalsglaeser, Material:Glas Hydrolytische Klasse III, braun, 500mL, DIN 55, Verschluesse fuer braune und weiBe Weithalsglaeser, schwarz, PP, LKD, DIN55	6pcs. Per Carton	3.P14802.500A
50mL white wide mouth glass jar, neck DIN 40, DIN 40 black PPN-caps, with white 2 mm PE-foam inlays	24pcs. Per Carton	3.P14803.0500
50mL amber wide mouth glass jar, neck DIN 40, DIN 40 black PPN-caps, with white 2 mm PE-foam inlays	24pcs. Per Carton	3.P14803.050A

Boston bottle

- suitable for storing liquid samples.
- suitable for storing light-sensitive samples.
- With a polypropylene cap and 0.015" PTFE liner.



Description	Packaging	Cat. No.
Boston Flauml;schchen, Material:Glas Hydrolytische Klasse III, braun, 30mL, neck 20-400, mit wei szlig;em Schraubverschluss, 0,010 "PTFE-Inlay, slash;33x78mm	24pcs. Per Carton	3.T12251.030A
Boston bottle, Material:Glas Hydrolytische Klasse III, braun, 60mL, neck 20-400, mit wei szlig;em Schraubverschluss, 0,010 "PTFE-Inlay, slash;39x94mm	24pcs. Per Carton	3.T12251.060A
Boston Flauml;schchen, Material:Glas Hydrolytische Klasse III, braun, 120mL, neck 22-400, mit wei szlig;em Schraubverschluss, 0,010 "PTFE-Inlay, slash;48x133mm	12pcs. Per Carton	3.T12251.120A
Boston Flaeschchen, Material:Glas Hydrolytische Klasse III, braun, 240mL, neck 28-400, mit weiBem Schraubverschluss, 0,010 "PTFE-Inlay, slash;Ø60x138mm	12pcs. Per Carton	3.T12251.240A
Boston Flauml;schchen, Material:Glas Hydrolytische Klasse III, braun, 480mL, neck 28-400, mit wei szlig;em Schraubverschluss, 0,010 "PTFE-Inlay, slash;75x169mm	6pcs. Per Carton	3.T12251.480A


Vials Rack



Description	Packaging	Cat. No.
Vails Rack for 12mm, 2mL vials, made from PP, white.	5 pcs. Per Carton	3.GSCM01.W005

Other colours products are on require.

Glass magnet



- Glass magnet is made from ploy urethane, which have good adsorptivity. So this product can adhere to the table surface and the samples without other assistant,
- This product is skidproof, flame resistant, antibacterial, easy to use, easy to transfer.
- This products solved the problem that sample vials are too small to operate, avoid losing or sloping when hand sampling in chromatogram laboratory.

Description	Packaging	Cat. No.
Glass magnet, blue, 160x90mm	1 pcs. Per PE bag	VHAP-1690-1
Glass magnet, blue, 160x90mm	5 pcs. Per Carton	VHAP-1690