

Crimping Tools and Decapping tools



- Coating out of the handle, ensure the tools can prevent hull corrosion;
- The head of the tools made from Special alloy materials ensure the long life;
- The handle part use frosted design ensure easy clench;
- Can be adjust the screw nut on the handle to ensure the best sealing performance;
- Crimping tools can use hexagon spanner to adjust the crimp neck height;
- Decapper easy to use, avoid the samples spilling and the glass parts broken;

Crimping tools

Description	Packaging	Cat. No.
Crimper for 8mm Aluminum caps	1 pc. Per box	3.009300.0800
Crimper for 11mm Aluminum caps	1 pc. Per box	3.009300.1100
Crimper for 13mm Aluminum caps	1 pc. Per box	3.009300.1300
Crimper for 20mm Aluminum caps	1 pc. Per box	3.009300.2000

Decapping tools

Description	Packaging	Cat. No.
Decapper for 8mm Aluminum caps	1 pc. Per box	3.009320.0800
Decapper for 11mm Aluminum caps	1 pc. Per box	3.009320.1100
Decapper for 13mm Aluminum caps	1 pc. Per box	3.009320.1300
Decapper for 20mm Aluminum caps	1 pc. Per box	3.009320.2000

Standard screw thread neck vials

These large volume vials is best choice for laboratory to storage the samples.

15-425 standard screw thread neck 8mL vials



15-425 screw thread neck vials

Description	Packaging	Cat. No.
15-425 Thread screw neck vial, 8ml, 61x16.6mm, clear glass, Borosilicate Type I Class A	100pcs. per Carton, 8 Cartons per Case	3.038015.0000
15-425 Thread screw neck vial, 8ml, 61x16.6mm, amber glass, Borosilicate Type I Class B	100pcs. per Carton, 8 Cartons per Case	3.038015.000A
15-425 Thread screw neck vial, 8ml, 61x16.6mm, clear glass, white marking spot and CNW LOGO, Borosilicate Type I Class A	100pcs. per Carton, 8 Cartons per Case	3.038015.00E0
15-425 Thread screw neck vial, 8ml, 61x16.6mm, amber glass, white marking spot and CNW LOGO, Borosilicate Type I Class B	100pcs. per Carton, 8 Cartons per Case	3.038015.00EA

15-425 screw thread neck preassembled caps

Description	Packaging	Cat. No.
Preassembled cap and septa for 15-425 thread screw, PP cap, black, centre hole, Clear silicone/White PTFE, 0.060" thick	100pcs. per PE bag	3.005350.1500
Preassembled cap and septa for 15-425 thread screw, PP cap, black, closed, Clear silicone/White PTFE, 0.060" thick	100pcs. per PE bag	3.005360.1500
Preassembled cap and septa for 15-425 thread screw, PP cap, black, centre hole, Red Butyl/Grey PTFE, 0.060" thick	100pcs. per PE bag	3.0F5350.1500
Preassembled cap and septa for 15-425 thread screw, PP cap, black, closed, Red Butyl/Grey PTFE, 0.060" thick	100pcs. per PE bag	3.0F5360.1500
Preassembled cap and septa for 15-425 thread screw, PP cap, white, closed, F217 line	100pcs. per PE bag	3.W05360.1500

15-425 screw thread neck caps

Description	Packaging	Cat. No.
15-425 Screw thread cap, made from PP, black, centre hole	100pcs. per PE bag	3.005310.1500
15-425 Screw thread cap, made from PP, black, closed	100pcs. per PE bag	3.005320.1500

